

Options & Specialty Items



MVS-1Y

MVS-1Y Flag Impeller Spindle

Use with the Small Sample Adapter to help keep sample materials in suspension



4B2

4B2 Spindle

Required for viscosity testing in accordance with ASTM D2983 (Low Temperature Viscosity Measurement of Automotive Fluid Lubricants)



ABZ

ABZ Spindle

Used for viscosity testing of thick film pastes. Short spindle length is suitable for immersion into shallow depth containers. Sensing length of spindle is less than 1-inch.

Custom Spindles

Custom spindles can be developed to meet your particular test requirements. Contact Brookfield or an authorized dealer for details.



SP1-UC-Y

YDX-1



SP1-UC-Y



EXTENSION LINK



SP-7Y QUICK CONNECT COUPLING

SP-4

SP-3

SP-5

Type D Extension Link

SPINDLE	TYPE
YDX-1	Male coupling nut; attaches to spindle
SP1-UC-Y	Female coupling nut; attaches to viscometer
SXV-XX	Extension link; see below for lengths

Type S Extension Link

SPINDLE
Available in LV-4 and RV/HA/HB spindle configurations

Extension

LINK	LENGTH	USED WITH
SXV-08	1"	UL Adapter
SXV-09	1.12"	Small Sample Adapter
SXV-24	3"	Thermosel
SXV-32	4"	Type D/S Extension
SXV-48	6"	Type D/S Extension
SXV-80	10"	Type D/S Extension
SXV-96	12"	Type D/S Extension

Note: Other lengths available; call for details

Quick Connect Coupling (SP-7Y)

PART	
SP-3	Coupling to viscometer/rheometer
SP-4	Coupling to spindle
SP-5	Sleeve (to hold together)



Labomat Essor
37 Bld Anatole France
F - 93287 Saint Denis Cedex
Tel.: +33 1 48 09 66 11
Fax: +33 1 48 09 98 65
E-mail : info@labomat.com
www.labomat.eu

Labomat Essor
Vlamingstraat 4
B - 8560 Wevelgem
Tel.: +32 56 43 28 13
Fax: +32 56 43 28 14
E-mail : info@labomat.com
www.labomat.eu



L.E Solutions
6 Imm B Résidence Ibn Batoua
Place Prince Sidi Mohammed
Belvédère
MA - 20300 Casablanca
Tel.: +212 52 22 41 714
Fax: +212 52 22 42 751
younesbaou@menara.ma
www.labomat.eu